

# PTHC125G



60" x 23.63" x 23.75" (H x W x D)

Single chamber passive heating blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, stainless steel cabinet, and lock

## Highlights:

Designed and assembled in the USA

No internal fans or moving parts ensure silent, safe, and efficient operation

Stainless steel construction inside and out for durable use and easy cleaning

## Product Features:

90 to 140°F temperature range

Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities

Passive heating

Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use

Digital thermostat

Electronic controls are conveniently located for easy temperature adjustment

Electronic readout

Large LED temperature display ensures your staff can easily monitor the interior temperature

Silent operation

No moving parts to ensure a quiet performance

Double pane tempered glass door

Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior

Professional handle

Stainless steel handle completes the professional look

Self-closing door



Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures

Stainless steel interior

Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning

## PTHC125G Specifications:

Over view	
Nominal Height	60.0" (152 cm)
Height of Cabinet	60.0" (152 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.75" (60 cm)
Depth with Handle	25.63" (65 cm)
Depth with door at 90°	46.0" (117 cm)
Capacity	10.57 cu.ft. (299 L)
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	150.0 lbs. (68 kg)
Parts & Labor Warranty	1 Year
UPC	761101066127
Warming Cabinet	
Door Swing	
Reversible	RHD
Adjustable Shelves	No
Shelf Type	Yes
Shelf Qty	Chrome
Thermostat Type	4
Interior Height	Digital
Interior Width	45.25" (115 cm)
Interior Depth	20.0" (51 cm)
Dolly	16.5" (42 cm)
Temperature Range	Optional
	90 to 140°F



Puretherm Warming Cabinet