

PTHC155G



71.5" x 23.63" x 23.75" (H x W x D)

Single chamber passive heating blanket warmer with a self-closing glass door, stainless steel interior, digital thermostat, black cabinet, and lock

Highlights:

Designed and assembled in the USA

No internal fans or moving parts ensure silent, safe, and efficient operation

Stainless steel interior ensures durable use and easy cleaning

Product Features:

90 to 140°F temperature range

Passive heating

Digital thermostat Electronic

readout Silent operation Double

pane tempered glass door

Professional handle

Self-closing door

Stainless steel interior

Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities

Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use

Electronic controls are conveniently located for easy temperature adjustment

Large LED temperature display ensures your staff can easily monitor the interior temperature

No moving parts to ensure a quiet performance

Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior



Stainless steel handle completes the professional look

Door will automatically close when left ajar, preventing items from accidentally being left exposed to ambient temperatures

Interior walls are constructed from stainless steel to ensure lasting durability and easy cleaning

PTHC155G Specifications:

Over view	
Nominal Height	72.0" (183 cm)
Height of Cabinet	71.5" (182 cm)
Nominal Width	24.0" (61 cm)
Width	23.63" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.75" (60 cm)
Depth with Handle	25.63" (65 cm)
Depth with door at 90°	46.75" (119 cm)
Capacity	12.17 cu.ft. (345 L)
Door	Glass
Cabinet	Black
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	170.0 lbs. (77 kg)
Parts & Labor Warranty	1 Year
UPC	761101066134
Warming Cabinet	
Door Swing	
Reversible	RHD
Adjustable Shelves	No
Shelf Type	Yes
Shelf Qty	Chrome
Thermostat Type	6
Interior Height	Digital
Interior Width	57.0" (145 cm)
Interior Depth	20.0" (51 cm)
Dolly	16.5" (42 cm)
Temperature Range	Optional
	90 to 140°F



Puretherm Warming Cabinet