

# PTHC65GUC



34.25" x 23.5" x 23.5" (H x W x D)

Undercounter single chamber passive heating blanket warmer with glass door, stainless steel interior and exterior, digital thermostat, and lock

## Highlights:

Designed and assembled in the USA

Temperature range: 90° to 140°F

No internal fans or moving parts ensure silent, safe, and efficient operation

## Product Features:

Undercounter Warming Cabinet

24" wide blanket warmer designed and constructed in the USA

90 to 140°F temperature range

Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities

Passive heating

Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use

Built-In capable

Front-breathing design lets you make the best use of space by installing your appliance under the counter

Digital thermostat

Electronic controls are conveniently located externally for easy temperature adjustment

Electronic readout

Large LED temperature display ensures your staff can easily monitor the interior temperature

Silent operation

No moving parts to ensure a quiet performance

Double pane tempered glass door

Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior

Stainless steel handle completes the professional look

Professional handle

## PTHC65GUC Specifications:

Over view	
Nominal Height	35.0" (89 cm)
Height of Cabinet	34.25" (87 cm)
Nominal Width	24.0" (61 cm)
Width	23.5" (60 cm)
Nominal Depth	24.0" (61 cm)
Depth	23.5" (60 cm)
Depth with Handle	25.38" (64 cm)
Depth with door at 90°	45.75" (116 cm)
Capacity	4.57 cu.ft. (129 L)
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	90.0 lbs. (41 kg)
Parts & Labor Warranty	1 Year
UPC	761101108056
Warming Cabinet	
Door Swing	
Reversible	RHD
Adjustable Shelves	Factory Reversible
Shelf Type	Yes
Shelf Qty	Chrome
Thermostat Type	3
Interior Height	Digital
Interior Width	26.0" (66 cm)
Interior Depth	19.38" (49 cm)
Level Legs	17.38" (44 cm)
Dolly	4
Temperature Range	Optional
	90 to 140°F

