

HORIZON[®] AIRGLIDE[®] RAIL HYDRAULIC STRETCHER



462 HORIZON® SERIES AIRGLIDE® RAIL HYDRAULIC STRETCHER

The Horizon Hydraulic AirGlide™ Series Stretcher is a multi-use unit to support a patient during induction, treatment, and recovery functions while avoiding unnecessary patient transfers to other equipment.



MOBILITY:

- Easy to maneuver: built-in central steering and braking systems with pedals conveniently located at all four corners
- Easy to position hydraulic foot pedals on both sides of the stretcher base allow for simple, hands-free height adjustment and Trendelenburg position
- Easy-to-operate casters and brakes to help reduce safety risks to healthcare personnel

PATIENT COMFORT:

- Continuous attention to the patient is possible with our unique Zone-of-Care™ design that allows all functions to be performed at the patient's side
- Unit allows for full length patient protection with side rails that are higher and longer allowing safe patient transfers with zero transfer gaps
- Backrest can be easily and consistently set at inclines from 0° - 90° with our cable-free, pneumatic-assisted fowler
- Knee-flex capability from 0° - 70°
- The base is equipped with an IV insert, storage space and a built in oxygen tank holder

462 Product Info

Specifications: 462 AIRGLIDE

	inch	cm [†]
Overall Length:	83	210
Overall Height: <i>range - minimum</i>	22	56
Overall Height: <i>range - maximum</i>	34.75	88
side rail Length:	75	191
side rail Height:	15	38
Tente® <i>easy rolling</i> Casters:	8	20
	degree of movement™	
Cable Free Fowler: <i>pneumatic-assist</i>	0 - 90	
Backrest: <i>movement</i>	0 - 90	
Knee Flex: <i>movement</i>	0 - 70	
Trendelenburg: <i>standard</i>	0 - 16	
Trendelenburg: <i>reverse</i>	0 - 16	
	lb	kg [†]
Shipping Weight:	370	168

Features:

Rail Type:	AirGlide™
IV Wells:	4 - 1/corner
Accessory Wells:	4 - 1/corner
Brake/Steer Pedals:	4 - 1/corner
Hydraulic Controls:	Foot Pedal - 1/side



426APAST PICTURED cover: with back tilted up, one rail up top; supine position, rails down far right; rails up



	462APAST		462CPAST		462DPAST	
	inch	cm [†]	inch	cm [†]	inch	cm [†]
Overall Width:	31.5	80	31.5	80	35	89
patient surface Width:	28	71	28	71	32	81
patient surface Length:	76	193	76	193	77	196
	feature		feature		feature	
Knee Gatch:	No		Yes		Yes	
	lb	kg [†]	lb	kg [†]	lb	kg [†]
Weight Capacity: <small>EVENLY DISTRIBUTED</small>	500	227	500	227	800	363

■ Standard
 ■ Optional
 — Does Not Apply
[†] same specs for both AC/110 volt and the DC/220 volt versions
[†] standard rounding method used to show metric values as whole numbers

GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
 Hausted is a registered trademark of GF Health Products, Inc.