

# Recovery Stretcher



## HS001

### Features:

- Locking mechanism with center-controlled locking system
- Steel Frame with aluminum alloy side-rails
- Built-in IV pole
- Trendelenburg: 20°
- Reverse Trendelenburg: 20°
- Backrest max angle: 75°
- Footrest max angle: 40°
- 330 lbs. capacity
- Mattress included; only available in black colour



\*Comes with black mattress only

### Specifications:

- Height: 25" - 38"
- Length: 84.5"
- Width: 31.6"

### Shipping Dimensions

- Height: 36"
- Length: 87"
- Width: 34"
- Weight: 170 kg.