



M-series MDS3030B06+MD30-EMG1

M-SERIES TALL EMERGENCY CART with MD30-EMG1 Emergency Package

Tall, Standard Width, Six Drawer Emergency Cart w/Breakaway Lock

Designed and built to meet the needs of today's healthcare market. Ergonomic side handles, comfort drawer pulls, expanded lock options and a standard plastic top make things easier for your staff and help support your infection control initiatives. M-Series carts offer you choices.

MD30-EMG1 Accessories Include:

Color shown is Red. More than 20 additional colors are available.



Drop shelf (MD-DRPSHLF)



Oxygen tank holder (O2HLDR)



IV pole (MD-IV2)



Cardiac board and brackets (MD-CARDBRD)



(2 Sets) Standard 3" Drawer Dividers (MD30-DIV3-S)



(2) Utility hooks (UHOOKS2)



Features:

- Sturdy, uni-body steel construction
- Choice of lock options
- Pack of 100 break-away seals
- 30" of vertical drawer space
- Full extension drawers slides for easy access
- Cart is delivered assembled
- Removable plastic top
- Latex free
- Fully customizable drawer configurations
- Industry best 5" quiet rolling casters - one steering caster and two braking casters
- Heat sealed powder coat steel break-away lock bar
- Ergonomic side push handles can be mounted on either side
- Body has pre-threaded accessory mounting holes
- Industry only fully curved protective bumper
- Durable scratch resistant heat sealed powder coat finish



Engineered & Manufactured in the U.S.

Drawer Configuration:

- Three 3" drawers (7.6 cm), two 6" drawers (15.2 cm) and one 9" drawer (22.7 cm)

Warranty: 12-Year Manufacturer's Warranty, 1-Year Limited Warranty on any bought-out parts

All dimensions are available in the "MA-Dims-30wide-5caster.pdf" document. Shipping weights can be found at the online product listing for this model, or feel free to request them from your Harloff Team.



View similar products online at

<https://www.harloff.com/product-category/emergency-crash-carts/>

Toll Free: (800) 433-4064 • Ph: (719) 637-0300
medicalsales@harloff.com • www.harloff.com